



# STEREO ZOOM MICROSCOPE

Stereo Zoom Microscopes produces a positive image. They are adopted with horizontal knob for zoom control; they have long working distance, clear and smooth image quality, wide field, convenience operation etc. They can be used for culture and education, scientific research, farming and forestry, electronics and precision machine industries widely.

- ▲ **VALUE**
- ▲ **VERSATILITY**
- ▲ **PERFORMANCE**

## MODEL : BSZ-41T



Model	BSZ-41	BSZ-41T
Viewing Head	Binocular	Trinocular
Eyepieces	WF 10x (20mm)	
Eyepiece tubes	45° and rotated 360°	
Objectives	1x to 4x	
Zoom Ratio	1:4	
Total Magnification	10x to 40x (Zoom)	
Field of View (mm)	23 to 5.5	
Working Distance (mm)	91	
Working Stage	Ø 95mm Frosted Glass	
Illumination	Incident: 12V, 10W	
	Transmitted : 12V, 10W	