



STEREO ZOOM MICROSCOPE

BSZ-45 Series Zoom Stereo Microscope is a kind of Binocular stereo microscope, which can magnify micro objects continually and show stereo up-right images. It provides clear high-contrast image, wide field and long working distance. It can be used for observation studies in medical and health, farming and forestry, as well as public security departments, schools and scientific research institutes, and is also used for inspection, assembling and repair of tiny spare parts in electronics and precision machine industries.

- ▶ **VALUE**
- ▶ **VERSATILITY**
- ▶ **PERFORMANCE**



MODEL : BSZ-45

TECHNICAL SPECIFICATIONS :

- > Magnification range : 7x to 45x
- > Objectives range if zoom magnification: 0.7X~4.5x
- > Eyepieces field is wide and clear : 10x (field size: Φ 20mm)
- > The binocular eyepiece tubes is inclined 45° and can be rotated 360°
- > The adjustment range of interpupilar distance: 53~75mm
- > The adjustment diopter of eyepiece tubes: -5~+5
- > Anti-mildew device is installed in the binocular to extend the life of instrument.
- > Transformer is out of the base of instrument
- > The mode of illumination selection and brightness adjusting are very easy, has three mode: transmitted light, reflected light and mixed light illumination
- > Different CCD camera adapter can be selected.
- > .5x LENS
- > LED RING ILLUMINATOR
- > ACROBAT STAND