

STEREO ZOOM MICROSCOPE

BSZ-5B



Description:

Stereo Zoom Microscope is a kind of binocular or trinocular stereo microscope, which can magnify micro objects continually and show up-right images.

Features:

- High quality, high function and excellent cost performance.
- Best choose with ergonomic design, sharp image display and wide viewing field.
- Providing a wide viewing field of $\Phi 20\text{mm}$ with the zoom ratio of 1:4.5.

BANBROS ENGINEERING PVT. LTD.

49/39, Site IV Industrial Area, Sahibabad,
Ghaziabad - 201010, Uttar Pradesh
Phone: +91 - 0120 - 4556371-76
Email: info@banbrosindia.com, sales@banbrosindia.com

Technical Specification:

Specifications		BSZ-5B	
Drawing Tube	Binocular head, inclined at 45°	●	
	Binocular head, inclined at 60°	○	
	Trinocular head, inclined at 45°	○	
	Digital head, inclined at 45°	○	
Diopter Adjustment Eyepiece	WF10X/Φ 20	●	
	WF15X/Φ16	○	
	WF20X/Φ12	○	
	WF30X/ Φ8	○	
	EW10X Eyepiece with cross hair	○	
Objective	Zoom objective	0.7X-3X	○
		1X-4.5X	●
	Auxiliary objective	2x WD: 30mm	○
		0.7X WD:120mm	○
		0.63X WD:137mm	○
	0.5X WD:180mm	○	
Zoom Ratio	1:4.3	○	
	1:4.5	●	
Working Distance	97mm	●	
Illumination	LED Reflection illumination	●	
	LED Transmission illumination	●	

Note:

- ● Standard Outfit
- ○ Optional

BANBROS ENGINEERING PVT. LTD.

49/39, Site IV Industrial Area, Sahibabad,
Ghaziabad - 201010, Uttar Pradesh
Phone: +91 - 0120 - 4556371-76
Email: info@banbrosindia.com, sales@banbrosindia.com