



**BANBROS ENGINEERING PVT. LTD.**  
The Precision Measurement

# METALLOGRAPHIC SPECIMEN GRINDING & POLISHING MACHINE

- *VALUE*
- *VERSATILITY*
- *PERFORMANCE*



**MBP-2BP**



#### DESCRIPTION:

Model **MBP-2B** metallographic specimen grinding & polishing machine is equipped with double disc and features stepless speed changing grinding and polishing.

The rotating speed of grinding and polishing disc is adjustable from 50 to 1000 rpm. It can realize the whole process of rough grinding, fine grinding, rough polishing and finishing polishing for specimen preparing.

This machine features easy operation and economical applying. It is the ideal specimen preparing instrument for using in factories, scientific research institutes and laboratories of colleges.

#### TECHNICAL SPECIFICATIONS:

Model No.	MBP-2BP
No. of Polishing disc	Double disc
Polishing head	With polishing head
Display of rotating speed	Digital display
Diameter of polishing disc	Ø 200mm
Diameter of grinding disc	Ø230mm
Polishing disc rotation	Variable
Polishing disc rotating speed	50-1000 RPM
Motor	550W 220V 50Hz
Dimensions	690 X 715 X 310mm
Net weight	45 Kg (approx.)

\* Image is for reference purpose only, due to continuous product development and upgradation the image and specification are subject to change.