



BANBROS ENGINEERING PVT. LTD.

The Precision Measurement People

PORTABLE METALLURGICAL MICROSCOPE



- *VALUE*
- *VERSATILITY*
- *PERFORMANCE*

PMM

**DESCRIPTION:**

This model is compact in size, light in weight & easy to carry around. This microscope permits examination on the desk & also on the post examination of an object in original size & shape in factories, laboratories & pipelines. The microscope unit is fitted with an incident and light through epi-illuminator, slot for dropping filter and Chargeable Battery Pack.

Two filters in mount (Green & Blue), & following optical combination in wooden box.

TECHNICAL SPECIFICATION:

- Eyepieces: 10x
- Objectives: 5x, 10x, 20 X & 40x
- Magnification: Upto 400X

*** Due to continuous product development, Image & specification can be upgrade.**