



STEREOZOOM MICROSCOPE



BSZ-608T

Features

- The The leading role of the same level microscopes, wide range observation with zoom range 0.6X-5X,
- zoom ratio 1:8.3 and visual field up to $\phi 76\text{mm}$.
- With high quality optical components, providing sharp and striking contrast image and no ghost image.
- Ease the tiredness with ergonomic design.
- h With Complete eyepieces, auxiliary objectives and other accessories, it is an ideal instrument for
- industry assembling, inspection and teaching field.

BANBROS ENGINEERING PVT. LTD.

49/39, Site IV Industrial Area, Sahibabad,
Ghaziabad - 201010, Uttar Pradesh
Phone: +91 - 0120 - 4556371-76
Email: info@banbrosindia.com, sales@banbrosindia.com

The Precision Measurement People

www.banbrosindia.com

TECHNICAL SPECIFICATIONS

Model	BSZ-608T
Viewing Tube	30° Trinocular viewing head, Inclined at 30°
Diopter Adjustment Eyepiece	WF10X/Ø23 Extra wide field eyepiece WF10X/Ø23
Zoom Objective	0.6X-5X
Magnification	6X to 50X
Zoom Ratio	1:8.3
Working Distance	115mm
Interpupillary Distance	55-75mm
Stage Insert	Glass plate, Diameter 125mm Opal & black acrylic plate, Diameter 125mm
Focusing Range	105mm
illumination	100V-240V/LED Incident Illumination 100V-240V/LED Transmitted Illumination
Video Adapter	1x C Mount

BANBROS ENGINEERING PVT. LTD.

49/39, Site IV Industrial Area, Sahibabad,
Ghaziabad - 201010, Uttar Pradesh
Phone: +91 - 0120 - 4556371-76
Email: info@banbrosindia.com, sales@banbrosindia.com

The Precision Measurement People

www.banbrosindia.com