



# METALLOGRAPHY SPECIMEN GRIDING & POLISHING MACHINE



## **BP-2G Series**

Model **BP-2G** is a cabinet type Metallography Specimen Polishing Machine. This machine features nice appearance design, big store cabinet and cooling system (which can clear up the heat produced during grinding so as to prevent the metallography structure from being damaged).

By equipping the motor with variable rotating speed, it can meet the different demands of different materials when polishing. It is the ideal specimen preparing instrument for using in factories, scientific research institutes and laboratories of colleges.

### TECHNICAL SPECIFICATIONS: -

Model No.	BP-2G	BP-2G(V)
No. of Polishing Disk	Double	Double
Diameter of polishing disk	200mm	200mm
Rotating Disk Type	Fixed	Variable
Polishing disk rotating speed	900 rpm	50-900 rpm
Motor	2x80W 220V 50Hz.	
Net weight	43 Kgs	

\* Due to continues product development, product image & Specifications are subject to change.

\* Customization can be done as per user requirement.