



STEREO ZOOM MICROSCOPE

DESCRIPTION: -

BSZ-5B Stereo Zoom Microscope is a kind of binocular or trinocular stereo microscope, which can magnify micro objects continually and show up-right images.

FEATURES: -

- High quality, high function and excellent cost performance.
- Best choose with ergonomic design, sharp image display and wide viewing field.
- Providing a wide viewing field of $\Phi 20\text{mm}$ with the zoom ratio of 1:4.5.



BSZ-5B

TECHNICAL SPECIFICATION: -

Specifications		BSZ-5B	
Drawing Tube	Binocular head, inclined at 45°	●	
	Binocular head, inclined at 60°	○	
	Trinocular head, inclined at 45°	○	
	Digital head, inclined at 45°	○	
Diopter Adjustment Eyepiece	WF10X/ Φ 20	●	
	WF15X/ Φ 16	○	
	WF20X/ Φ 12	○	
	WF30X/ Φ 8	○	
	EW10X Eyepiece with cross hair	○	
Objective	Zoom objective	0.7X-3X	○
		1X-4.5X	●
	Auxiliary objective	2x WD: 30mm	○
		0.7X WD:120mm	○
		0.63X WD:137mm	○
	0.5X WD:180mm	○	
Zoom Ratio	1:4.3	○	
	1:4.5	●	
Working Distance	97mm	●	
Illumination	LED Reflection illumination	●	
	LED Transmission illumination	●	

Note: ● Standard Outfit ○ Optional