



**BANBROS ENGINEERING PVT. LTD.**

The Precision Measurement People

# DIGITAL VIDEO LCD MICROSCOPE

- *VALUE*
- *VERSATILITY*
- *PERFORMANCE*



**BSZ-450 LCD**



**DESCRIPTION:**

**BSZ-450 LCD** Series are new generation Stereo zoom video microscope with microcomputer which drives the operation system Win CE5.0 and display system. It provides clear high-contrast image, wide field and long working distance.

It can be used for observation studies in medical and health, farming and forestry, as well as public security departments, schools and scientific research institutes, and is also used for inspection, assembling and repair of tiny spare parts in electronics and precision machine industries.

**FEATURES:**

- The Image can be easily viewed by many people on LCD monitor at the same time.
- Improve the working condition also can soften the tiredness caused by a long-time observation.
- Integrated structure, easily operation, Small size and can save space.
- The Magnification can be zoomed by non-step.
- There is “+” for purpose of positioning and measurement.
- Easily equipment, it can be work properly by put it on stable place.

**TECHNICAL SPECIFICATIONS:**

Model No.	BSZ-450 LCD
LCD Screen	8" inch
Reduce Lens	0.5X
Objective	0.65X to 4.5X
CCD Resolution	1/3" 580 Line
Incident Light	12V/10W Halogen Lamp
Transmitted Light	12V/10W Halogen Lamp

\* Due to continuous product development, Image & specification can be upgrade.