



BANBROS ENGINEERING PVT. LTD.

The Precision Measurement People

VICKER HARDNESS TESTER



- *VALUE*
- *VERSATILITY*
- *PERFORMANCE*

BTV-250/250-S

**FEATURES:**

- Table top sturdy model
- Front projection screen with micrometer with optical system list count of 0.0005 mm
- Motorised loading & unloading of loads
- Load selection by indexible load selection knob
- Very wide testing range from very soft metals to very hard metals
- Models with minimum 1 kg load are available
- The illumination light goes off during the application of load which reduces the operator's fatigue due to continuous light illumination as well as increases bulb life
- State of art load damping system for a very smooth & jerk free loading & unloading of applied load
- Machine confirms to is: 1754-2002 for vickers & is: 2281-2005 / bs240 for brinell.
- Four model are available :
- BTV-50 with loads 5 to 50 kgf
- BTV-250 with loads 5 to 250 kgf
- BTV-50 (spl.) With loads 1 to 50 kgf
- BTV250 (spl.) With loads 1 to 250 kgf

TECHNICAL SPECIFICATIONS:

Model No.	BTV-250	BTV-250-S
Load (Kgf)	5, 10, 20 ,30, 50 for Vickers 62.5, 187.5, 250 for Brinell	1, 5, 10, 20, 30, 50 for Vickers 62.5, 187.5, 250 for Brinell
Maximum Test height (mm)	200 mm	
Depth of throat (mm)	150 mm	
Weight Approx. (Kgf)	180	
Dimension L x W x H mm	660 x 250 x 850 mm	
Connected load	0.1 KVA, 1 Phase AC, 220V, 50 Hz	

STANDARD ACCESSORIES:

- Flat Table of 100mm Dia: 1 no.
- Vee Table for Specimen- 6 to 80 mm Dia: 1 no.
- Vickers Diamond Indenter: 1 no.
- Ø2.5, 5mm Ball Indenter: 1 no. each.
- Ø2.5, 5mm Steel Balls : 1 no. each.
- Vickers test block: 1 no.
- Brinell test block: 1 no.
- Electrical Cord: 1 no.
- Instruction Manual: 1 no.