



METALLOGRAPHY SPECIMEN GRINDING & POLISHING MACHINE

Model **MBP Series** Metallographic Specimen grinding & polishing machine is equipped with single disc and features stepless speed changing grinding and polishing.

The rotating speed of grinding and polishing disc is adjustable from 50 to 1000 rpm. It can realize the whole process of rough grinding, fine grinding, rough polishing and finishing polishing for specimen preparing.

This machine features easy operation and economical applying. It is the ideal specimen preparing instrument for using in factories, scientific research institutes and laboratories of colleges.



MBP-1B



MBP-2B

TECHNICAL SPECIFICATIONS: -

Model No.	MBP-1B	MBP-2B
No. of Polishing Disk	Single	Double
Display of Rotating Speed	Digital Display	Digital Display
Diameter of polishing disk	Ø 200mm	Ø 200mm
Diameter of grinding disk	Ø230mm	Ø230mm
Polishing disk Rotation	Variable	Variable
Polishing disk rotating speed	50-1000 RPM	50-1000 RPM
Motor	550W 220V 50Hz	
Dimensions	65 x 37 x 31 cm	690 X 715 X 310mm
Net weight	30 Kgs	45 Kg.